



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **γ-Secretase Inhibitor III**

Cat. number: **VD-PI-0014**


CAS Number: **Not available.**

Lot number: **VD13-0777105**

Personal protective equipment

■ Lab coat.

■ Safety glasses.



Product Specification

Formula:

C₂₀H₃₀N₂O₄

MW:

362.5

Purity:

≥95%.

Appearance:

Off-white solid.

Solubility:

DMSO.

Storage:

-20°C. **Hygroscopic.**

Synonyms:

N-Benzylloxycarbonyl-Leu-leucinal; Z-LL-CHO

Product Description:

A cell-permeable, reversible inhibitor of gamma-secretase.

Signature _____