



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **L-C-Propargylglycine**

Cat. number: **VD-Ant-0341**

CAS Number: **23235-01-0**


Lot number: **VD13-06631100**

Personal protective equipment

■ Lab coat.



■ Safety glasses.



Product Specification

Formula: C₅H₇NO₂

MW: 113.12

Purity: >=99%

Appearance: White powder.

Solubility: DMSO.

Storage: +2°C to +8°C. **Protect from light.**

Synonyms: H-Pra-OH

Product Description:

This unusual amino acid is produced by microorganisms, it has been detected in the poisonous mushroom *Amanita abrupta*. Precursor for tritiated norleucine which may be introduced synthetically into peptides.

Product Reference:

- 1.P.Scannell et al., J. Antibiot., 24, 239 (1971).
- 2.O.Leukart et al., Helv. Chim. Acta, 59, 2181 (1976).

Signature _____