



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **Hexanoyl Glycine**

Cat. number: **VD-NP-0095**

CAS Number: **24003-67-6**


Lot number: **VD12-0047908**

Personal protective equipment

■ Lab coat.



■ Safety glasses.



Product Specification

Formula:	C ₈ H ₁₅ NO ₃
MW:	173.21
Purity:	>=95% by HPLC
Appearance:	White solid
Solubility:	DMSO
Storage:	-20°. Protect from light.
Synonyms:	NSC 224460; Caproylglycine

Product Description:

An acylaminoacid derivative which is used as an anti-psychotic drug.

Product Reference:

1. Willard, H., et al.: J. Clin. Invest., 65, 690 (1980).

Signature _____