



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com


CERTIFICATE OF ANALYSIS

Product name: **Potassium D-tartrate monobasic**

Cat. number: **VD-NP-00041** CAS Number: **57341-16-9** Lot number: **VD13-00005721**

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₄ H ₅ KO ₆
MW:	188.18
Purity:	>=98%
Appearance:	White powder or crystals.
Solubility:	50mg/ml, H ₂ O.
Storage:	+2°C to +8°C. Protect from light.
Synonyms:	Potassium 3-carboxy-2,3-dihydroxypropanoate

Product Description:

Not available.

Product Reference:

1.Dawson, R.M.C. Data for Biochemical Research 3rd ed., Oxford , (1986), 421.

Signature _____