



MANUFACTURED BY VDM BIOCHEMICALS, INC.  
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA  
PHONE NUMBER: 1-440-786-9400  
FAX NUMBER: 1-440-786-9410  
EMAIL: [SALES@VDMBIO.COM](mailto:SALES@VDMBIO.COM)

[www.vdmbio.com](http://www.vdmbio.com)


## CERTIFICATE OF ANALYSIS

Product name: **L-Cycloserine**

Cat. number: **VD-OEIA-0157** CAS Number: **339-72-0; 4834-58-6** Lot number: **VD12-009988**

### Personal protective equipment

- Lab coat.
- Safety glasses.



### Product Specification

Formula:	C <sub>3</sub> H <sub>6</sub> N <sub>2</sub> O <sub>2</sub>
MW:	102.1
Purity:	>=98% by HPLC.
Appearance:	White to off-white powder.
Solubility:	Water.
Storage:	-20°. <b>Hygroscopic. Unstable in neutral or acidic solutions.</b>
Synonyms:	S(-)-4-Amino-3-isoxazolidinone
NMR:	<b>Consistent with structure.</b>

### Product Description:

Excitatory amino acid. Inhibitor of ketosphinganine synthase.

### Quality controle department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE