



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS


Product name: **Ac-LETD-AFC**

Cat. number: **VD-CRP-0217** CAS Number: **210345-02-1** Lot number: **VD-CP60039**

www.vdmbio.com

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₃₁ H ₃₈ F ₃ N ₅ O ₁₂
MW:	729.6
Purity:	>=97%
Appearance:	White solid.
Solubility:	DMSO.
Storage:	-20°C. Protect from light.
Synonyms:	Ac-Leu-Glu-Thr-Asp-AFC
NMR:	Consistent with structure.

Product Description

Specific fluorogenic substrate for caspase-8.

Quality controle department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE