



MANUFACTURED BY VDM BIOCHEMICALS, INC.  
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA  
PHONE NUMBER: 1-440-786-9400  
FAX NUMBER: 1-440-786-9410  
EMAIL: SALES@VDMBIO.COM

[www.vdmbio.com](http://www.vdmbio.com)


## CERTIFICATE OF ANALYSIS

Product name: **Ac-WEHD-AFC**

Cat. number: **VD-CRP-0189** CAS Number: **Not available.** Lot number: **VDAA-000053** [www.vdmbio.com](http://www.vdmbio.com)

### Personal protective equipment

- Lab coat.
- Safety glasses.



### Product Specification

<b>Formula:</b>	C <sub>38</sub> H <sub>37</sub> F <sub>3</sub> N <sub>8</sub> O <sub>11</sub>
<b>MW:</b>	838.74
<b>Purity:</b>	≥95%
<b>Appearance:</b>	Yellow solid.
<b>Solubility:</b>	DMSO or H <sub>2</sub> O.
<b>Storage:</b>	+2°C to +8°C. <b>Protect from light.</b>
<b>Synonyms:</b>	Ac-Trp-Glu-His-Asp-AFC
<b>NMR:</b>	<b>Consistent with structure.</b>

### Product Description

Fluorogenic substrate for caspase 1 and 5.

### Quality control department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE