



MANUFACTURED BY VDM BIOCHEMICALS, INC.  
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA  
PHONE NUMBER: 1-440-786-9400  
FAX NUMBER: 1-440-786-9410  
EMAIL: [SALES@VDMBIO.COM](mailto:SALES@VDMBIO.COM)

[www.vdmbio.com](http://www.vdmbio.com)

## CERTIFICATE OF ANALYSIS

Product name: **Endoxifen**


Cat. number: **VD-ACR-0014**

CAS Number: **112093-28-4**

Lot number: **VD15-117008**

### Personal protective equipment

- Lab coat.
- Safety glasses.



### Product Specification

<b>Formula:</b>	C <sub>25</sub> H <sub>27</sub> NO <sub>2</sub>
<b>MW:</b>	373.49
<b>Purity:</b>	≥98%.
<b>Appearance:</b>	A crystalline solid.
<b>Solubility:</b>	DMSO to 100 mM.
<b>Storage:</b>	+2°C to +8°C.
<b>Synonyms:</b>	Not available.
<b>LC-MS:</b>	Conforms to structure.

### Product Description

Metabolite of tamoxifen.

### Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE