



MANUFACTURED BY VDM BIOCHEMICALS, INC.
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **Y-27632 (free base)**


Cat. number: **VD-IS-0134**

CAS Number: **146986-50-7**

Lot number: **VD16-03320**

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₁₄ H ₂₁ N ₃ O
MW:	247.34
Purity:	≥98%.
Appearance:	Solid.
Solubility:	Water.
Storage:	Store at -20°C.
Synonyms:	Rock Inhibitor

Product Description

A highly potent, cell-permeable, reversible, and selective inhibitor of Rho-associated protein kinases.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE