



MANUFACTURED BY VDM BIOCHEMICALS, INC.
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **Necrostatin-1, Inactive Control**

Cat. number: **VD-AIH-0125**


CAS Number: **64419-92-7**

Lot number: **VD14-08034**

Personal protective equipment

■ Lab coat.

■ Safety glasses.



Product Specification

Formula: C₁₂H₁₁N₃O₃S

MW: 245.3

Purity: ≥95% by HPLC.

Appearance: Light yellow solid

Solubility: DMSO or MeOH.

Storage: Stock solutions are stable for up to 3 months at -20°C.

Synonyms: 5-(Indol-3-ylmethyl)-2-thiohydantoin

Product Description

A cell-permeable N-demethylated thiohydantoin analog of Nec-1 that is devoid of anti-necroptotic properties and serves as a suitable inactive control.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE