



MANUFACTURED BY VDM BIOCHEMICALS, INC.
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com


CERTIFICATE OF ANALYSIS

Product name: **MU-PHE-HOMOPHE-Fluoromethylketone**

Cat. number: **VD-AIH-0285** CAS Number: **Not available.** Lot number: **VD15-005135**

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₂₅ H ₃₀ FN ₃ O ₄
MW:	455.52
Purity:	≥98%.
Appearance:	White solid.
Solubility:	DMSO.
Storage:	Store at +4°C.
Synonyms:	Mu-Phe-Homophe-FMK; Mu-F-HomoPhe-FMK

Product Description

Cathepsin B and L inhibitor peptide.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE