



MANUFACTURED BY VDM BIOCHEMICALS, INC.  
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA  
PHONE NUMBER: 1-440-786-9400  
FAX NUMBER: 1-440-786-9410  
EMAIL: SALES@VDMBIO.COM

[www.vdmbio.com](http://www.vdmbio.com)

## CERTIFICATE OF ANALYSIS

Product name: **L-GLY-7-Amino-4-Methylcoumarin**


Cat. number: **VD-CRP-0207**

CAS Number: **N/A**

Lot number: **VD14-049001**

### Personal protective equipment

- Lab coat.
- Safety glasses.



### Product Specification

<b>Formula:</b>	C12H12N2O3
<b>MW:</b>	232.23
<b>Purity:</b>	>=98%.
<b>Appearance:</b>	Off-white solid.
<b>Solubility:</b>	DMSO or H2O.
<b>Storage:</b>	Store at -20°C.
<b>Synonyms:</b>	L-Gly-AMC; L-G-AMC

### Product Description

N/A.

### Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE