



MANUFACTURED BY VDM BIOCHEMICALS, INC.
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **AY 9944 dihydrochloride**


Cat. number: **VD-OEIA-0058**

CAS Number: **366-93-8**

Lot number: **VD16-0809341**

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₂₂ H ₂₈ Cl ₂ N ₂ · 2HCl
MW:	464.30
Purity:	≥99%
Appearance:	White solid.
Solubility:	50 mM in Water.
Storage:	Desiccate at RT.
Synonyms:	Trans-N,N-bis[2-Chlorophenylmethyl]-1,4-cyclohexanedime thanamine dihydrochloride

Product Description

Inhibitor of hedgehog (Hh) signaling.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE