



MANUFACTURED BY VDM BIOCHEMICALS, INC.
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **MNS**

Cat. number: **VD-IS-0079**

CAS Number: **1485-00-3**


Lot number: **VD14-22209**

Personal protective equipment

■ Lab coat.



■ Safety glasses.



Product Specification

Formula:	C ₉ H ₇ NO ₄
MW:	193.16
Purity:	≥99%.
Appearance:	Yellow solid.
Solubility:	100 mM in DMSO and to 5 mM in Ethanol.
Storage:	Store at +4°C.
Synonyms:	3,4-Methylenedioxy-b-nitrostyrene

Product Description

Selective inhibitor of Src and Syk tyrosine kinases.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE