



MANUFACTURED BY VDM BIOCHEMICALS, INC.

5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA

PHONE NUMBER: 1-440-786-9400

FAX NUMBER: 1-440-786-9410

EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **VEGFR2 Kinase Inhibitor II**

Cat. number: **VD-GFC-0122**


CAS Number: **288144-20-7**

Lot number: **VD15-45798**

Personal protective equipment

■ Lab coat.

■ Safety glasses.



Product Specification

Formula: C₁₇H₁₅BrN₂O

MW: 343.2

Purity: ≥95%.

Appearance: Brown solid.

Solubility: DMSO.

Storage: -20°C.

Synonyms: Vascular Endothelial Growth Factor Receptor 2 Kinase Inhibitor II

Product Description

A membrane permeant indolin-2-one class of receptor tyrosine kinase (RTK) inhibitor. Does not inhibit EGF-R kinase activity.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE