



MANUFACTURED BY VDM BIOCHEMICALS, INC.
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com


CERTIFICATE OF ANALYSIS

Product name: **Z-VAL-DL-ASP-Fluoromethylketone**

Cat. number: **VD-CRP-0212** CAS Number: **582316-00-5** Lot number: **VD14-10501**

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula: C₁₈H₂₃FN₂O₆

MW: 382.38

Purity: ≥98%.

Appearance: White powder.

Solubility: DMSO.

Storage: Store at -20°C.

Synonyms: Z-VD-FMK

Product Description

An irreversible general caspase inhibitor.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE