



MANUFACTURED BY VDM BIOCHEMICALS, INC.

5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA

PHONE NUMBER: 1-440-786-9400

FAX NUMBER: 1-440-786-9410

EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **L-Leucine 7-amido-4-methylcoumarin hydrate**

Cat. number: **VD-OEIA-0029**


CAS Number: **66447-31-2**

Lot number: **VD16-4700099**

Personal protective equipment

■ Lab coat.

■ Safety glasses.



Product Specification

Formula: C₁₆H₂₀N₂O₃

MW: 288.37

Purity: ≥95%.

Appearance: Solid.

Solubility: DMSO.

Storage: +2°C to +8°C.

Synonyms: H-LEU-AMC

Product Description

This product is sensitive fluorogenic substrate for leucine aminopeptidase and aminopeptidase M.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE