



MANUFACTURED BY VDM BIOCHEMICALS, INC.
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **D-VLK-AMC**

Cat. number: **VD-PI-0180**


CAS Number: **Not available.**

Lot number: **VD14-00020**

Personal protective equipment

■ Lab coat.

■ Safety glasses.



Product Specification

Formula:	C27H41N5O5
MW:	515.64 (Not including 2TFA salt).
Purity:	>=98%.
Appearance:	White solid.
Solubility:	DMSO.
Storage:	RT.
Synonyms:	D-VAL-LEU-LYS-7-amino-4-methylcoumarin

Product Description

Plasmin, tissue plasminogen activator substrate.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE