



MANUFACTURED BY VDM BIOCHEMICALS, INC.
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com


CERTIFICATE OF ANALYSIS

Product name: **DL-Glyceraldehyde 3-phosphate diethyl acetal barium salt**

Cat. number: **VD-OEIA-0209** CAS Number: **93965-35-6** Lot number: **VD15-22240**

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C7H15BaO7P
MW:	379.49
Purity:	Not available.
Appearance:	White powder.
Solubility:	Clear to slightly hazy colorless solution at 50 MG plus 1 ML of 1 N HCL.
Storage:	Store at -20°C.
Synonyms:	3,3-Diethoxy-1,2-propanediol 1-phosphate barium salt
Proton NMR spectrum:	CONFORMS TO STRUCTURE.
Water by Karl Fischer:	< OR = 2.0%

Product Description

An inhibitor of aspartate aminotransferase isozymes.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE