



MANUFACTURED BY VDM BIOCHEMICALS, INC.
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **Ac-VEID-FMK**

Cat. number: **VD-CRP-0224**


CAS Number: **Not available**

Lot number: **VD14-0077**

Personal protective equipment

■ Lab coat.

■ Safety glasses.



Product Specification

Formula:	C ₂₅ H ₄₁ FN ₄ O ₉
MW:	560.6
Purity:	Single spot by TLC.
Appearance:	Lyophilized.
Solubility:	DMSO.
Storage:	Store at +4°C.
Synonyms:	Ac-V-E(OMe)-I-D(OMe)-FMK

Product Description

A cell-permeable, irreversible inhibitor of caspase-6,-8,-10.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE