



MANUFACTURED BY VDM BIOCHEMICALS, INC.
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **FITC-YVAD(OMe)-FMK**


Cat. number: **VD-CRP-0225**

CAS Number: **Not available**

Lot number: **VD15-0099000**

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula: C₄₄H₄₄FN₅O₁₂S

MW: 886.0

Purity: >=98%

Appearance: Orange solid.

Solubility: DMSO.

Storage: 2-8°C.

Synonyms: Not available.

MS: Consistent with structure.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE