



MANUFACTURED BY **VDM BIOCHEMICALS, INC.**
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **Rho/SRF Pathway Inhibitor, CCG-1423**

Cat. number: **VD-GFC-0126**


CAS Number: **285986-88-1**

Lot number: **VD14-1110**

Personal protective equipment

■ Lab coat.

■ Safety glasses.



Product Specification

Formula: C₁₈H₁₃ClF₆N₂O₃

MW: 454.8

Purity: ≥97% by HPLC.

Appearance: White solid.

Solubility: DMSO (10mg/ml).

Storage: Store at +4°C.

Synonyms: CCG-1423

LC-MS: Confirms the structure

Product Description

A cell-permeable benzamide compound that acts as an effective inhibitor against RhoA- and RhoC-mediated cellular activities by targeting signaling events.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE