



MANUFACTURED BY VDM BIOCHEMICALS, INC.
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **CHC**


Cat. number: **VD-OEIA-0149**

CAS Number: **28166-41-8**

Lot number: **VD16-27000**

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

| | |
|--------------------|---|
| Formula: | C10H7NO3 |
| MW: | 189.17 |
| Purity: | >=98%. |
| Appearance: | Yellow solid. |
| Solubility: | 100 mM in DMSO and to 50 mM in Ethanol. |
| Storage: | +2°C to +8°C. Protect from light. |
| Synonyms: | (2E)-2-cyano-3-(4-hydroxyphenyl)prop-2-enoic acid |

Product Description

Monocarboxylic acid transport (MCT) inhibitor.

Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE