



MANUFACTURED BY VDM BIOCHEMICALS, INC.  
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA

PHONE NUMBER: 1-440-786-9400  
FAX NUMBER: 1-440-786-9410  
EMAIL: [SALES@VDMBIO.COM](mailto:SALES@VDMBIO.COM)

[www.vdmbio.com](http://www.vdmbio.com)

## CERTIFICATE OF ANALYSIS

Product name: **Z-WEHD-FMK**


Cat. number: **VD-CRP-0213**

CAS Number: **210345-00-9 (TFA salt)**.

Lot number: **VD15-05005**

### Personal protective equipment

- Lab coat.
- Safety glasses.



### Product Specification

<b>Formula:</b>	C37H42FN7O10
<b>MW:</b>	763.78
<b>Purity:</b>	>=95%.
<b>Appearance:</b>	White lyophilized solid.
<b>Solubility:</b>	DMSO (15.28 mg/ml).
<b>Storage:</b>	-20°C.
<b>Synonyms:</b>	Benzyloxycarbonyl-Trp-Glu(OMe)-His-Asp(OMe)-fluoromethylketone

### Product Description

Caspase-5 inhibitor. Inhibits ECyd-induced rRNA fragmentation in mouse mammary tumor and human fibrosarcoma cells.

### Quality control department:

This material is chemically synthesized. Made in USA

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE