



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com


CERTIFICATE OF ANALYSIS

Product name: **Caspase-9 Inhibitor III**

Cat. number: **VD-CRP-0181** CAS Number: **403848-57-7** Lot number: **VD16-1615601** www.vdmbio.com

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₂₄ H ₃₅ ClN ₆ O ₉
MW:	587.0
Purity:	>=95% by HPLC
Appearance:	White solid.
Solubility:	DMSO or Water
Storage:	-20°C
Synonyms:	Ac-LEHD-CMK
NMR:	Consistent with structure.

Product Description

Supplied as a trifluoroacetate salt. A potent, irreversible inhibitor of caspase-9. Reported to reduce myocardial infarct size during reperfusion (~70 nM).

Quality control department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE