



MANUFACTURED BY VDM BIOCHEMICALS.LLC  
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA  
PHONE NUMBER: 1-440-786-9400  
FAX NUMBER: 1-440-786-9410  
EMAIL: SALES@VDMBIO.COM

[www.vdmbio.com](http://www.vdmbio.com)


## CERTIFICATE OF ANALYSIS

Product name: **(R)-2-Acetylthio-3-phenylpropionic acid**

Cat. number: **VD-OEIA-00091** CAS Number: **57359-76-9** Lot number: **VD12-4100909**

### Personal protective equipment

- Lab coat.
- Safety glasses.



### Product Specification

<b>Formula:</b>	C <sub>11</sub> H <sub>12</sub> O <sub>3</sub> S
<b>MW:</b>	224.28
<b>Purity:</b>	>=95% by HPLC
<b>Appearance:</b>	Orange oil
<b>Solubility:</b>	DMSO
<b>Storage:</b>	-4°C. <b>Protect from light.</b>
<b>Synonyms:</b>	D-2-Thioacetyl-3-phenylpropionic Acid

### Product Description:

A IMP-1 metallo-β-lactamase inhibitor.

### Product Reference:

1.Greenlee, M.L., et al.: Boorg. Med. Chem. Lett., 9, 17, 2549 (1999).

*Signature* \_\_\_\_\_