



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **2-Guanidinoethylmercaptosuccinic acid**


Cat. number: **VD-OEIA-1104**

CAS Number: **77482-44-1**

Lot number: **VD776654000**

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₇ H ₁₃ N ₃ O ₄ S
MW:	235.2
Purity:	>=99% by TLC
Appearance:	White solid
Solubility:	H ₂ O.
Storage:	+2°C to +8°C. Protect from light.
Synonyms:	GEMSA
NMR:	Consistent with structure.

Product Description

A potent enkephalin convertase inhibitor.

Quality controle department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE