



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com


CERTIFICATE OF ANALYSIS

Product name: **2-Amino-2-norbornanecarboxylic acid**

Cat. number: **VD-OEIA-1106** CAS Number: **20448-79-7** Lot number: **VD00008766** www.vdmbio.com

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

| | |
|-------------|------------------------------------------------|
| Formula: | C ₈ H ₁₃ NO ₂ |
| MW: | 155.19 |
| Purity: | >=98% |
| Appearance: | White to off-white powder. |
| Solubility: | 50mg/ml in 1N NH ₄ OH. |
| Storage: | +2°C to +8°C. Protect from light. |
| Synonyms: | BCH |
| NMR: | Consistent with structure. |

Product Description

Inhibitor of amino acid transport across cell membranes.

Quality controle department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE