



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com


CERTIFICATE OF ANALYSIS

Product name: **N,N-Dimethyl Sunitinib**

Cat. number: **VD-AM-0130** CAS Number: **326914-17-4** Lot number: **VD9088274** www.vdmbio.com

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₂₀ H ₂₃ FN ₄ O ₂
MW:	370.42
Purity:	>=95%
Appearance:	Solid.
Solubility:	Water.
Storage:	RT. Keep cool and dry.
Synonyms:	Not available.
NMR:	Consistent with structure.

Product Description

N,N-Dimethyl Sunitinib is a Sunitinib intermediate that is an angiogenesis inhibitor and antitumor prodrug.

Quality controle department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE