



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS


Product name: **Sirtinol**

Cat. number: **VD-OEIA-0081** CAS Number: **410536-97-9** Lot number: **VDAA-065433** www.vdmbio.com

Personal protective equipment

■ Lab coat.

■ Safety glasses.



Product Specification

Formula:	C ₂₆ H ₂₂ N ₂ O ₂
MW:	394.5
Purity:	>=97%
Appearance:	Off-white to yellow crystals.
Solubility:	DMSO.
Storage:	-20°C. Protect from light.
Synonyms:	Sir Two Inhibitor Naphthol
NMR:	Consistent with structure.

Product Description

Sirtinol is a cell permeable specific inhibitor of sirtuin NAD-dependant deacetylases.

Quality controle department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE