



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **TPEN**

Cat. number: **VD-AI-0269**


CAS Number: **16858-02-9**

Lot number: **VD-000122**

www.vdmbio.com

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₂₆ H ₂₈ N ₆
MW:	424.6
Purity:	>=97%
Appearance:	Solid.
Solubility:	Dimethyl formamide, DMSO, 100% Ethanol or Methanol.
Storage:	+2°C to +8°C . Protect from light.
Synonyms:	N,N,N',N'-Tetrakis-(2-pyridylmethyl)ethylenediamine
NMR:	Consistent with structure.

Product Description

A cell permeable and extremely high-affinity chelator of heavy metals, yet with low affinity for Mg²⁺ and Ca²⁺.

Quality controle department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE