



MANUFACTURED BY VDM BIOCHEMICALS.LLC  
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA  
PHONE NUMBER: 1-440-786-9400  
FAX NUMBER: 1-440-786-9410  
EMAIL: SALES@VDMBIO.COM

[www.vdmbio.com](http://www.vdmbio.com)


## CERTIFICATE OF ANALYSIS

Product name: **Letrozole**

Cat. number: **VD-Ant-0256** CAS Number: **112809-51-5** Lot number: **VDAA-007230** [www.vdmbio.com](http://www.vdmbio.com)

### Personal protective equipment

- Lab coat.
- Safety glasses.



### Product Specification

<b>Formula:</b>	C <sub>17</sub> H <sub>11</sub> N <sub>5</sub>
<b>MW:</b>	285.31
<b>Purity:</b>	>=98%
<b>Appearance:</b>	White or off-white crystalline powder.
<b>Solubility:</b>	Dichloromethane.
<b>Storage:</b>	+2°C to +8°C. <b>Protect from light.</b>
<b>Synonyms:</b>	CGS-20267
<b>NMR:</b>	<b>Consistent with structure.</b>

### Product Description

A nonsteroidal aromatase inhibitor structurally related to Fadrozole, that has been introduced for the adjuvant treatment of hormonally-responsive breast cancer.

### Quality controle department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE