



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS


Product name: **Z-LLL-AMC**

Cat. number: **VD-OEIA-0159** CAS Number: **Not available.** Lot number: **VD-CP009138**

www.vdmbio.com

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₃₆ H ₄₈ N ₄ O ₇
MW:	648.8
Purity:	>=98%.
Appearance:	White lyophilized powder.
Solubility:	DMSO.
Storage:	-20°C. Protect from light.
Synonyms:	Z-Leu-Leu-Leu-AMC
NMR:	Consistent with structure.

Product Description

Fluorogenic proteasome substrate.

Quality controle department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE