



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com


CERTIFICATE OF ANALYSIS

Product name: **Z-IETD-FMK**

Cat. number: **VD-NP-0071** CAS Number: **Not available.** Lot number: **VDAA-000041** www.vdmbio.com

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₃₀ H ₄₃ FN ₄ O ₁₁
MW:	654.68
Purity:	>=98%
Appearance:	Yellow solid
Solubility:	DMSO
Storage:	-20°C. Protect from light.
Synonyms:	Caspase-8 inhibitor
NMR:	Consistent with structure.

Product Description

Irreversible, cell permeable inhibitor of caspase-8 and granzyme B. Inhibits also caspase-10. The inhibitor is designed as a methyl ester to facilitate cell permeability.

Quality control department: _____

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE