



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com


CERTIFICATE OF ANALYSIS

Product name: **Ambenonium dichloride**

Cat. number: **VD-OEIA-1110** CAS Number: **52022-31-8** Lot number: **VD-11044317** www.vdmbio.com

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₂₈ H ₄₂ Cl ₄ N ₄ O ₂
MW:	608.47
Purity:	>=98%
Appearance:	White solid.
Solubility:	H ₂ O.
Storage:	-20°C. Protect from light.
Synonyms:	Not Available.
NMR:	Consistent with structure.

Product Description

Not available.

Quality controle department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE