



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com



CERTIFICATE OF ANALYSIS

Product name: **PDK1/Akt/Flt Dual Pathway Inhibitor**

Cat. number: **VD-AIH-0288** CAS Number: **331253-86-2** Lot number: **VD1000997** www.vdmbio.com

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₁₀ H ₄ N ₆ O
MW:	224.2
Purity:	>=98% (sum of two isomers).
Appearance:	Yellow solid.
Solubility:	DMSO
Storage:	+2°C to +8°C. Protect from light.
Synonyms:	Akt Inhibitor XXI, PDK1 Inhibitor I
NMR:	Consistent with structure.

Product Description

A cell-permeable compound that selectively induces apoptosis in AML (Acute Myelogenous Leukemia) with little effect on normal CD34+ AML progenitor cells.

Quality control department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE