



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS


Product name: **Ac-VEID-AMC**

Cat. number: **VD-CRP-0216** CAS Number: **219137-97-0** Lot number: **VD-CP60038**

www.vdmbio.com

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₃₂ H ₄₃ N ₅ O ₁₁
MW:	673.7
Purity:	>=95% by HPLC.
Appearance:	White to off-white powder.
Solubility:	DMSO; dilute with distilled water or buffer.
Storage:	-20°C. Protect from light.
Synonyms:	Ac-Val-Glu-Ile-Asp-AMC; Caspase-6 substrate (fluorogenic).
NMR:	Consistent with structure.

Product Description

Fluorogenic substrate for caspase-6 and related cysteine proteases.

Quality control department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE