



MANUFACTURED BY VDM BIOCHEMICALS.LLC  
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400  
FAX NUMBER: 1-440-786-9410  
EMAIL: SALES@VDMBIO.COM

[www.vdmbio.com](http://www.vdmbio.com)

## CERTIFICATE OF ANALYSIS

Product name: **γ-Secretase Inhibitor II**

Cat. number: **VD-PI-0181**


CAS Number: **Not available.**

Lot number: **VD13-0777100**

### Personal protective equipment

■ Lab coat.

■ Safety glasses.



### Product Specification

**Formula:** C<sub>33</sub>H<sub>57</sub>F<sub>2</sub>N<sub>5</sub>O<sub>9</sub>

**MW:** 705.8

**Purity:** >=90%.

**Appearance:** White solid.

**Solubility:** DMSO.

**Storage:** -20°C. **Protect from light.**

**Synonyms:** Mw167

### Product Description:

A cell-permeable, reversible and selective peptidomimetic inhibitor of γ-secretase

*Signature* \_\_\_\_\_