



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **Ac-IETD-CHO**


Cat. number: **VD-CRP-0214**

CAS Number: **191338-86-0**

Lot number: **VD13-0000983**

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₂₁ H ₃₄ N ₄ O ₁₀
MW:	502.5
Purity:	>=98%
Appearance:	Off-white solid.
Solubility:	DMSO.
Storage:	RT. Keep cool and dry.
Synonyms:	Ac-Ile-Glu-Thr-Asp-CHO

Product Description: _____
Caspase activator.

Product Reference: _____
1.Not available.

Signature _____