



MANUFACTURED BY VDM BIOCHEMICALS.LLC  
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400  
FAX NUMBER: 1-440-786-9410  
EMAIL: SALES@VDMBIO.COM

[www.vdmbio.com](http://www.vdmbio.com)


## CERTIFICATE OF ANALYSIS

Product name: **EBPC**

Cat. number: **VD-OEIA-0110** CAS Number: **57056-57-2** Lot number: **VD13-04447180**

### Personal protective equipment

- Lab coat.
- Safety glasses.



### Product Specification

<b>Formula:</b>	C <sub>14</sub> H <sub>15</sub> NO <sub>4</sub>
<b>MW:</b>	261.28
<b>Purity:</b>	>=99%
<b>Appearance:</b>	White solid.
<b>Solubility:</b>	DMSO to 100 mM.
<b>Storage:</b>	+2°C to +4°C. <b>Protect from light.</b>
<b>Synonyms:</b>	Ethyl 1-benzyl-3-hydroxy-2(5H)-oxopyrrole-4-carboxylate

### Product Description:

Potent aldose reductase inhibitor.

### Product Reference:

1.Mylari et al (1991).J.Med.Chem. 34 1011.

*Signature* \_\_\_\_\_