



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com


CERTIFICATE OF ANALYSIS

Product name: **Z-LEED-FMK**

Cat. number: **VD-CRP-0215** CAS Number: **Not available.** Lot number: **VD13-00005775**

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₃₂ H ₄₅ FN ₄ O ₁₂
MW:	696.7
Purity:	>=95%
Appearance:	Solid.
Solubility:	DMSO.
Storage:	-20°C. Protect from light.
Synonyms:	Z-Leu-Glu(OMe)-Glu(OMe)-Asp(OMe)-fluoromethyl ketone

Product Description: _____
Cell permeable Caspase-13 inhibitor.

Product Reference: _____
Not available.

Signature _____