



MANUFACTURED BY VDM BIOCHEMICALS.LLC  
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400  
FAX NUMBER: 1-440-786-9410  
EMAIL: SALES@VDMBIO.COM

[www.vdmbio.com](http://www.vdmbio.com)

## CERTIFICATE OF ANALYSIS

Product name: **Boc-Lys(Boc)-OH**


Cat. number: **VD-AA-0197**

CAS Number: **2483-46-7**

Lot number: **VD12-05772314**

### Personal protective equipment

- Lab coat.
- Safety glasses.



### Product Specification

<b>Formula:</b>	C <sub>16</sub> H <sub>30</sub> N <sub>2</sub> O <sub>6</sub>
<b>MW:</b>	346.4
<b>Purity:</b>	>=99%
<b>Appearance:</b>	White powder
<b>Solubility:</b>	DMSO.
<b>Storage:</b>	+2°C to +8°C. <b>Keep cool and dry.</b>
<b>Synonyms:</b>	(S)-2,6-Bis-tert-butoxycarbonylaminohexanoic acid

### Product Description:

Not available.

### Product Reference:

Not available.

*Signature* \_\_\_\_\_