



MANUFACTURED BY VDM BIOCHEMICALS.LLC

5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400

FAX NUMBER: 1-440-786-9410

EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **N-(4-Carboxy-3-hydroxyphenyl)maleimide**

Cat. number: **VD-OEIA-1103**


CAS Number: **19232-43-0**

Lot number: **VD13-1199954**

Personal protective equipment

■ Lab coat.

■ Safety glasses.



Product Specification

Formula:

C₁₁H₇NO₅

MW:

233.182

Purity:

≥99%

Appearance:

Solid.

Solubility:

DMSO.

Storage:

RT. **Keep cool and dry.**

Synonyms:

4-Maleimidylsalicylic acid

Product Description:

Reported to be an irreversible inhibitor of LDH isoenzyme.

Product Reference:

1.Carvajal, G., et al., National Cancer Institute Monograph, 27: 111, (1967).

Signature _____