



MANUFACTURED BY VDM BIOCHEMICALS.LLC  
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400  
FAX NUMBER: 1-440-786-9410  
EMAIL: SALES@VDMBIO.COM

[www.vdmbio.com](http://www.vdmbio.com)

## CERTIFICATE OF ANALYSIS

Product name: **Scriptaid**


Cat. number: **VD-IS-0087**

CAS Number: **287383-59-9**

Lot number: **VD011472**

### Personal protective equipment

- Lab coat.
- Safety glasses.



### Product Specification

**Formula:** C<sub>18</sub>H<sub>18</sub>N<sub>2</sub>O<sub>4</sub>

**MW:** 326.35

**Purity:** >=99% by TLC

**Appearance:** Yellow solid.

**Solubility:** DMSO to 100 mM.

**Storage:** -20°C. **Protect from light.**

**Synonyms:** Not available.

### Product Description:

Novel histone deacetylase inhibitor. Produces >100-fold increase in histone acetylation and displays relatively low toxicity.

### Product Reference:

1.Su et al (2000) A novel histone deacetylase inhibitor identified by high-throughput transcriptional screening of a compound library. Cancer Res. 60 3137.

*Signature* \_\_\_\_\_