



MANUFACTURED BY VDM BIOCHEMICALS.LLC
5386 MAJESTIC PARKWAY, UNIT-9,BEDFORD HEIGHTS, OH 44146,USA

PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **H-Asp(Gly-OH)-OH**


Cat. number: **VD-Ant-0333**

CAS Number: **3790-52-1**

Lot number: **VD12-0009718**

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₆ H ₁₀ N ₂ O ₅
MW:	190.16
Purity:	>=95%
Appearance:	White powder.
Solubility:	DMSO.
Storage:	-20°C. Protect from light.
Synonyms:	H-β-Asp-Gly-OH

Product Description:

Was found to be a substance unique to caecal contents of germ-free and antibiotic-treated mice. Was also detected in rat kidney, isolated from uremic serum and found in human urine.

Product Reference:

- 1.T.Abiko et al., Biochem. Biophys. Res. Commun., 82, 707 (1978).
- 2.G.W.Welling and G.Groen, Biochem. J., 175, 807 (1978).
- 3.T.Tanaka et al., J. Biochem. (Tokyo), 84, 1147 (1978).

Signature _____