



MANUFACTURED BY VDM BIOCHEMICALS, INC.
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **Autophagy Inhibitor, 3-MA**


Cat. number: **VD-OEIA-10020**

CAS Number: **5142-23-4**

Lot number: **VD76609103**

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₆ H ₇ N ₅
MW:	149.15
Purity:	>=99%
Appearance:	White solid
Solubility:	DMSO
Storage:	+2°C to +8°C. Protect from light.
Synonyms:	3-methyladenine
NMR:	Consistent with structure.

Product Description

A synthetic intermediate and a cell-permeable autophagic sequestration blocker that protects cerebellar granule cells from apoptosis post serum/potassium deprivation.

Quality control department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE