



MANUFACTURED BY VDM BIOCHEMICALS, INC.
5386 MAJESTIC PARKWAY, UNIT-9, BEDFORD HEIGHTS, OH 44146, USA
PHONE NUMBER: 1-440-786-9400
FAX NUMBER: 1-440-786-9410
EMAIL: SALES@VDMBIO.COM

www.vdmbio.com

CERTIFICATE OF ANALYSIS

Product name: **Mitomycin C**

Cat. number: **VD-AI-0272**


CAS Number: **50-07-7**

Lot number: **VD117124**

www.vdmbio.com

Personal protective equipment

- Lab coat.
- Safety glasses.



Product Specification

Formula:	C ₁₅ H ₁₈ N ₄ O ₅
MW:	334.33
Purity:	>=98%.
Appearance:	Blue-violet solid.
Solubility:	5 mM in H ₂ O and to 100 mM in DMSO.
Storage:	+2°C to +8°C. Protect from light.
Synonyms:	Ametycine; MMC
NMR:	Consistent with structure.

Product Description

Antibiotic and antitumor agent.

Quality controle department:

This material is chemically synthesized. **Made in USA**

TO BE USED FOR RESEARCH PURPOSES ONLY. NOT FOR DRUG OR CLINICAL USE IN HUMANS OR FOOD ADDITIVE USE